

# Documentation HeartRatePlugin ver. 14

Content:

- [Important Notice.....2](#)
- [Description.....2](#)
- [Requirements.....3](#)
- [Preferences.....4](#)
- [HeartRatePluginEvent.....6](#)
- [Model.....7](#)
- [Methods.....8](#)
- [Properties.....8](#)
- [How to use the HeartRatePlugin.....9](#)
- [Examples.....10](#)
- [Known issues.....11](#)

## **Important Notice:**

This version differs in structure basically to the previous versions!

## **Description**


This plugin scans, connects and reads data from a heart rate sensor(hrs), based on events.  
You can connect to multiple hrs at the same time

# Requirements


- Bluetooth and GPS have to be turned on
- The heart rate sensor(hrs) has to use the standard heart rate protocol defined by SIG:  
Heart Rate Service: 0x180D  
Heart Rate Measurement Characteristics: 0x2A37

How to check, if your device supports the standard heart rate protocol:

- Install nRF Connect app
  - connect to your hrs
- Heart Rate should be displayed




- Tap on Heart Rate
- Heart Rate Measurement should be displayed




## Preferences

Android:

- Check out Unity/Preferences
- In External Tools tab under Android please be sure to have all checkboxes activated



- Please check in Unity Hub, if all modules for Android Build Support are added.



✓	✓	Android Build Support	328.8 MB	1.0 GB
	✓	Android SDK & NDK Tools	1.0 GB	3.1 GB
	✓	OpenJDK	165.0 MB	72.7 MB

- Check, if the correct platform for the plugin is set:
  - select Assets>HeartRatePlugin>Scripts>Bridge>HeartRateMulti
  - In the Inspector under Select platforms for plugin, Android has to be selected (nothing else)

## MacOS:

- Check, if the correct platform for the plugin is set:
  - select Assets>HeartRatePlugin>Scripts>Bridge>MacOSHeartRatePlugin.bundle
  - In the Inspector under Select platforms for plugin, Editor and Standalone have to be selected (nothing else)

## iOS:

- Check, if the correct platform for the plugin is set:
  - select Assets>HeartRatePlugin>Scripts>Bridge>libiOSHeartRatePlugin
  - In the Inspector under Select platforms for plugin, iOS has to be selected (nothing else)

## HeartRatePlugin.Event:

### Properties:

- **MacID:**
  - Identifier(iOS) or MAC-Adress(Android) of the sensor, who raised the Event
- **Info:**
  - Information depending on the EventType
- **Type:**
  - **SYSTEM\_NOT\_SCANNING**
    - raised, if scan is stopped. Contains reason in Info
  - **SYSTEM\_SCANNING**
    - raised, if scan is started
  - **NEW\_SENSOR**
    - raised, if a new hrs was found upon scanning
  - **REMOVE\_UNCONNECTED**
    - raised, if a sensor is removed from HeartRateSensor-Model
  - **CONNECTING**
    - raised, if a sensor is connecting
  - **CONNECTED**
    - raised, if a sensor is connected
  - **DISCONNECTED**
    - raised, if a sensor got disconnected
  - **NOTIFICATION\_CONTROLPOINT**
    - raised, if Heart Rate Control Point is detected
  - **NOTIFICATION\_BODYSENSORLOCATION**
    - raised, if Body Sensor Location is detected
  - **NOTIFICATION\_MEASUREMENT**
    - raised, if a sensor send measurement data

## Model

The HeartRateSensor-model contains informations and data of the hrs

### Fields:

- **MacId:**
  - Identifier(iOS) or MAC-Adress(Android)
- **Name:**
  - Name the sensor provides
- **Rssi:**
  - Distance from hrs to MobileDevice
- **IsConnecting:**
  - System started a connection to the hrs
- **IsConnected:**
  - Connection was successfully established
- **PulseRate:**
  - Actual pulse rate the hrs send
- **SCStatus:**
  - Sensor contact status of the hrs:
    - NOT\_SUPPORTED
    - NOT\_SUPPORTED\_1
    - NO\_CONTACT
    - CONTACT
- **EnergyExpended:**
  - actual amount of expended energy. [Null](#) if not present
- **RrInterval:**
  - Array of last heart rate variability. [Null](#) if not present
- **SensorLocation:**
  - Body sensor location of the hrs. Mandatory.
    - OTHER
    - CHEST
    - WRIST
    - FINGER
    - HAND
    - EARLOBE
    - FOOT
- **HR\_ControlPoint:**
  - Control point of the Alert Notification server. [Null](#) if not supported

All discovered sensors are stored in the Dictionary [HeartRateSensor.Sensors\[key\]](#).  
Key: MacId.

## Methods

- **public void StartScan():**  
Checks if mobile device is capable to scan for bluetooth devices and starts the scan process.  
Sets `IsScanning` to `true`.  
At Start, all unconnected hrs are removed from `HeartRateSensor.Sensors`  
`EventType.SYSTEM_SCANNING` is raised.
- **public void StopScan()**  
Stops the scan process.  
Sets `IsScanning` to `false`.  
`EventType.SYSTEM_NOT_SCANNING` is raised.
- **public void Connect(string MacId)**  
Connects to the given `MacId`.  
Sets `HeartRateSensor.Sensors[MacId].IsConnecting` to `true`.  
`EventType.CONNECTING` is raised
- **public void Disconnect(string MacId, bool All)**  
Disconnects the given `MacId`  
`All`: all hrs connected to the mobile device will be disconnected

## Properties

- **IsInitialized**  
Returns `true` if the plugin is completely initialized
- **IsScanning**  
Returns `true` if system is scanning



## How to use the HeartRatePlugin

Please take a look at the example scenes.

1. Create an empty object in your scene
  2. It is mandatory to name it exactly HeartRatePlugin and attach HeartRatePlugin.cs to it (Assets/HeartRatePlugin/Scripts/HeartRatePlugin.cs)
- Setup for your scene script:
3. Create reference to HeartRatePlugin:  
[SerializeField]  
private HeartRatePlugin heartRatePlugin;
  4. Drag HeartRatePlugin-object to the property in the Inspector of your scene script
  5. Create method:  
void OnHeartRateEvent(object sender, HeartRatePlugin.EventArgs e) with switch(e.Type) as shown in SuperSimpleExample.cs
  6. In your Start() or OnEnable() (or where you want to), attach HeartRatePlugin.Event to this method:  
HeartRatePlugin.Event += OnHeartRateEvent;
  7. Start scanning with:  
heartRatePlugin.StartScan();
  8. React to the arriving Events:  
E.g. ff you get EventType.NEW\_SENSOR you can connect to it via:  
heartRatePlugin.Connect(e.MacId);
  9. Settings of the sensor arrive at:
    - EventType.NOTIFICATION\_CONTROLPOINT:  
check if Heart Rate Control Point is supported at  
HeartRateSensor.Sensors[e.MacId].HR\_ControlPoint.  
O or Null means NOT\_SUPPORTED
    - EventType.NOTIFICATION\_BODYSENSORLOCATION:  
check what sensor location the hrs is attached to at  
HeartRateSensor.Sensors[e.MacId].SensorLocation.
  10. Data arrives at EventType.NOTIFICATION\_MEASUREMENT:  
Using e.MacID you can access all data at HeartRateSensor.Sensors[e.MacId].Foo (→ Fields)

## Examples

Two Example scenes are included.

- **MultiConnection:**

This Example shows, what this plugin is capable of. Some visualization styles are included:



select visualization style



e.g. Bumping Heart



e.g. Ecg

- **SuperSimpleExample:**

Pure scan and connect, data is logged

Setup for SuperSimpleExample:

1. Create empty object
2. Name it as HeartRatePlugin
3. Attach HeartRatePlugin.cs to it (Assets/HeartRatePlugin/Scripts/HeartRatePlugin.cs)
4. Attach SuperSimpleExample.cs
5. Drag Object HeartRatePlugin to Script SuperSimpleExample
6. Take a look at the logs

## Known issues

- If you run unity on a windows machine, you can build for iOS and Andoid, but you can't play in editor.
- Please be sure, you choose the correct platform.